# FRAUNHOFER INSTITUTE FOR SOLAR ENERGY SYSTEMS ISE DIVISION »HYDROGEN TECHNOLOGIES«

**Ex-situ characterisation tools for materials and components in PEM water electrolysis and fuel cells** 



Dr. A. Georg, U. Groos, Dr. M. Klingele, Dr. S. Metz, Dr. S. Ouardi, Dr. A. Schaadt, Dr. T. Smolinka, M. Zedda


Fraunhofer Institute for Solar Energy Systems ISE

www.h2-ise.com

www.ise.fraunhofer.com



#### Ex situ analysis of materials and components Overview




Capillary flow porometry (CFP) for through-plane gas permeability of porous media

Porosimetry for absolute and relative in-plane gas and liquid permeability



High temperature and Near ambient pressure X-ray photoelectron spectroscopy (HT NAP-XPS)



Electrical conductivity and interfacial contact resistance (ICR) and forcedisplacement measurements



Environmental Scanning Electron Microscopy (ESEM) and Energy Dispersive X-ray Spectroscopy (EDS)



X-ray computed tomography (Micro CT)




Inductively coupled plasma mass spectroscopy (ICP –MS) for element analysis of liquids



# **Capillary flow porometry for porous media**

Gas-liquid displacement porometry with pressure step/stability method:

- Determination of pore diameter:
  - Largest / mean / smallest
- Through-plane gas permeability
- Relationship capillary pressure vs saturation
- Inner (mean) contact angle
- Important parameters for PTL design





**3** Bromberger, K., Ghinaiya, J., Lickert, T., Fallisch, A., Smolinka, T.: Hydraulic ex situ through-plane characterization of porous transport layers in PEM water electrolysis cells, Int. J. Hydrogen Energy, Vol. 43, Issue 5, 2018, 2556-2569, https://doi.org/10.1016/j.ijhydene.2017.12.042.



### In plane permeability of porous media

- Absolute and relative in-plane gas and liquid permeability
- Test set up with a cell area of 25 cm<sup>2</sup>
- Method established for GDL / PTL etc. using Darcy or Darcy-Forchheimer law






Test set up for in-plane permeability measurements and detailed view on test cell



#### Electric conductivity and interfacial contact resistance

Understanding of internal cell resistance and individual contributions is crucial for stack design


- Measurement of in-plane and through-plane resistance/conductivity as function of compression
- Interfacial Contact Resistance (ICR) of passivation layers
- Thickness measurement by eddy current sensors (resolution 3 μm)
- Combination of clamping pressure and thickness measurement allows a force/displacement analysis





#### Through-plane resistance and interfacial contact resistance

- Measurement of through-plane resistance of different materials
- Interfacial contact resistance (ICR) between bipolar plate and GDL
- Bulk resistance of GDL
- Thickness measurement






6

#### Inplane & through-plane resistance





Fraunhofer



#### **Electrochemical characterization of bipolar plate materials**

CV of 316L at different temperatures

EI: 0.001mol/I H,SO, + 0.1 mg/I HF

E = -0.1 bis 1.5 V (Cycle 2)

Scanrate: 0.1 mV/s

active corrosion

H<sub>o</sub> production

-0.2 0.0 0.2 0.4

AE: 316L

RF RHF

Ar-Saturation

GE: Pt

25 -

20


15 17/Cm<sup>2</sup> 10

current density /

-5

-10

-15



Test cell made of PTFE, integrated heating, gas (Ar or O<sub>2</sub>), working electrode (WE) (sample), reference electrode (RHE), counter electrode (CE). Cyclic voltammogramm of stainless steel at different temperatures, electrolyte:  $0.001 \text{ M H}_2\text{SO}_4+0.1 \text{ mg HF}$ 

0,6

passivation

O<sub>2</sub> reduction

0,8

voltage / V

1,0

current density 25°C cycle 2


current density 60°C cycle 2

current density 80°C cycle 2

1,2 1,4 1,6

transpassivation


O<sub>production</sub>



Electrochemical measurements is combined with elemental analysis of the electrolyte with ICP-MS (here during potentiostatic test at 0.8 V for 4 days), SEM/EDX analysis, and contact resistance measurement.



#### Inductively Coupled Plasma – Mass Spectroscopy (ICP-MS)




Element analysis of liquids with ICP-MS. In Ar-plasma the molecules are destroyed and ionised. In the mass spectrometer the ions of every mass are counted. ICP-MS analysis of electrolyte after ex-situ corrosion test (24 h at 0.8 V) of metallic bipolar plates. (ISE-coating compared with Au plated and uncoated stainless steel substrate)

ICP-MS product water analysis of a fuel cell stack with metallic bipolar plates (extreme example)



9

## **Scanning Electron Microscopy and Computed X-ray Tomography**




Zielke, L., Fallisch, A., Paust, N., Zengerle, R., Thiele S.: Tomography based screening of flow field/current 10 collector combinations for PEM water electrolysis. RSC Adv (2014), 10.1039/c4ra12402b © Fraunhofer



#### **Fenton testing**


- Investigating the chemical stability of fuel cell or electrolysis membranes regarding radical attacs
- Membrane degradation: cations catalyze radical formation from H<sub>2</sub>O<sub>2</sub> and radicals attac the membrane polymer
- Fenton test: Insertion of membrane in 30% H<sub>2</sub>O<sub>2</sub> with metal salts (typically Fe<sup>2+</sup>) and heat at 80°C
- Measurement of Fluoride-release over time as measure for membrane decomposition
- F- measurement with ion sensitive electrode
- Use of chemical inert polymer bottles (safer and less contamination than glass)
- Testing with different temperatures, cations and/or membranes; variation of cation concentration





#### Near ambient pressure X-ray photoelectron spectroscopy (NAP-XPS)

- X-ray source: AlKα (1.4 keV)
- Spot size 300 µm Ø
- Pressure:
  10<sup>-8</sup> mbar up to 20 mbar
- Sample temperature:
  5 800 °C
- Sample size up to 120 mm Ø
- Plasma cleaning of the surface
- Fully automated vacuum and gas dosing system (N<sub>2</sub>, Ar, H<sub>2</sub>, O<sub>2</sub>)



Surface characterisation of an IrO<sub>2</sub> based anode for PEM water electrolysis as received, after short-term and long-term operation EnvironESCA: Electron spectroscopy for chemical analysis under environmental conditions at near ambient pressure of catalysts, liquids and liquid-solid interfaces.



12

#### **Thanks a lot for your kind attention!**



Fraunhofer Institute for Solar Energy Systems ISE

**Electrolysis** 

Dr. Tom Smolinka tom.smolinka@ise.fraunhofer.de Fuel Cell Systems Ulf Groos ulf.groos@ise.fraunhofer.de Thermochemical Systems

Dr. Achim Schaadt achim.schaadt@ise.fraunhofer.de

